


# Endoscope Video Camera HD

**visionflex**  
Equality healthcare


Impressive HD image quality packaged in an extraordinary European design.

Light, ergonomic and versatile – the Endoscope Camera HD delivers high-precision HD images for state-of-the-art endoscopy. Focusing on the aspects of ergonomics, hygiene and cutting-edge technology, an entirely new camera concept was developed, meeting all the industry requirements.

This endoscope also excels in terms of size and weight – an ultra-small camera head weighing a mere 70 grams. As a result, the camera is barely noticeable within the endoscopic system itself.



**Flawless design.  
Ideal handling.**



# Endoscope Video Camera HD

Impressive HD image quality packaged in an extraordinary European design.

Thanks to its superb combination of low weight and small footprint, the imaging device boasts excellent handling in any situation. Incorporating the latest pioneering technology available, M-CAM HD can be operated without special drivers and with no need for installation.

The camera is simply connected to a USB 2.0 port.

The M-CAM HD is controlled via two freely configurable buttons located on the camera head. An LED status display within the camera head provides details of the function selected. As an option, all functions can also be displayed directly on the live video screen.

## Available in two options:

M-CAM HD (without integrated light source) or M-CAM HD+ (with integrated LED light source)

## FEATURES & BENEFITS

- High-resolution sensor
- Waterproof aluminum housing
- Innovative user-interface
- HD endo coupler
- Haptic feedback
- M-CAM HD+: Integrated LED light source



## SPECIFICATIONS

<b>Weight</b>	70 grams (without TV-Adapter)
<b>Dimensions</b>	41 x 40 x 32 mm (W x H x D)
<b>Housing</b>	aluminium, stainless steel
<b>Resolution</b>	HD 720p
<b>Single Images</b>	2952 x 1944 px
<b>Control</b>	2 freely configurable buttons
<b>Device Status</b>	LED status display – On-screen overlay
<b>Connection</b>	USB 2.0 port or higher
<b>Operating System</b>	Windows – macOS
<b>Update</b>	Online update function
<b>Driver</b>	UVC standard technology
<b>Disinfection</b>	IP 68 – 24-hour immersion
<b>ON/OFF</b>	Motion detection

For more information or to book a demonstration, please contact Visionflex at [sales@visionflex.com](mailto:sales@visionflex.com)



+61 2 8914 4000



[visionflex.com](http://visionflex.com)



**visionflex**  
Equality healthcare