

Video Examination Glasses HD

visionflex
Equality healthcare

Capturing clear and precise high-quality medical images has never been easier on this handheld device.

Ideal for telehealth and education these extremely lightweight HD video glasses offer a unique level of diagnostic insight by viewing the patient from the medical practitioner's perspective.

These high-performance, HD-image glasses are ideal for team collaboration, hands-on demonstrations and training, and real-time teleconferencing. The camera's high-resolution and sensitivity makes them ideal for use in a wide variety of settings.



**High definition,
hands free imaging**



Video Examination Glasses HD

HD video glasses for hands free imaging with USB connection.

Simple to use with collapsible temples and ergonomic nose pads.

Perfect for virtual care consultations, connecting to a ProEX system, Mac or PC through 3 metre USB cable.

SPECIFICATIONS

Built-in lens	1.83mm/F2.0 - untinted
Image sensor	1/6" Digital CMOS Camera
Image area	2.7288 x 1.5498mm
Pixel size	1.4 x 1.4
Resolution	1280 x 720 @30fps
Sensitivity	553 mV/(Lux-sec)
Storage temperature	-30°C to 60°C
Working temperature	-1°C to 45°C
Interface	USB 2.0 UVC
Power source	DC5V (Micro USB OTG B Type)
Power current	60mA
Dimension	160 x 150mm

FEATURES & BENEFITS

- High performance full HD image glasses
- 3 metre USB 2.0 cable with standard USB-A connector
- Supplied in protective box



For more information or to book a demonstration, please contact Visionflex at sales@visionflex.com

+61 2 8914 4000

 visionflex.com



visionflex
Equality healthcare