

Non-Contact Infrared Forehead Thermometer with Bluetooth

Provides accurate touch free temperature measurements in just one second

This non-contact forehead thermometer is a multifunctional device that provides accurate temperature measurement in just one second. Suitable for temperature readings of the body and surfaces, at a distance of 5cm. It allows for the non-invasive and gentle measurement of body temperature without having to touch a patient.

**Non-invasive,
touch-free
measurement**



Non-Contact Infrared Forehead Thermometer with Bluetooth

visionflex
Equality healthcare


The thermometer uses infrared technology to measure body and surface temperatures quickly and accurately. Measures body temperatures between 32°C (89.6°F) and 43°C (109.4°F), with an accuracy of $\pm 0.2^{\circ}\text{C}$ for temperatures between 35°C and 42°C and $\pm 0.03^{\circ}\text{C}$ for temperatures below 35°C or above 42°C. The high/low indicator and dual colour backlit display assists with quick diagnosis. Built in memory stores up to 30 most recent readings. Bluetooth connectivity enables transmission of data between the thermometer to devices including the Visionflex ProEX Telehealth Hub, PC or tablet.

Features & benefits

- **Suitable for measuring temperatures of adults, children and surfaces**
- **One touch measurement – fast and accurate results in one second**
- **Memory stores up to 30 readings for easy recall**
- **Fever indicator with high temperature alarm and flash backlight**
- **High/low indicator and dual colour backlit display**
- **Temperature reading in Fahrenheit or Celsius**
- **Bluetooth connectivity**
- **60 second timeout period**

Specifications

Displayed temperature range (body): 32°C to 43°C (89.6°F to 109.4°F)
Displayed temperature range (surface): 0°C to 100°C (32°F to 212°F)
Measurement accuracy (body): $\pm 0.2^{\circ}\text{C}$ for the range of 35.0°C to 42.0°C / $\pm 0.3^{\circ}\text{C}$ for the range of $< 35.0^{\circ}\text{C}$ or $> 42.0^{\circ}\text{C}$
Measurement accuracy (surface): $\pm 1^{\circ}\text{C}$
Weight: 61.8g (without battery)
Dimensions: 155 x 40 x 39mm
Power source: 2 x 1.5V AAA alkaline batteries
Battery life: approx. 3,000 measurements
Operating temperature: 10°C to 40°C (50°F to 104°F)
Operating humidity: 30% to 85% relative humidity (non-condensing)
Operating Atmosphere Pressure: 700hPa to 1060hPa
Memory capacity: 30 measurements
Water resistance: IP22
Timeout period: 60 seconds

 For further information and a demonstration, contact **Visionflex** at info@visionflex.com.au

 1300 059 926 +61 2 8438 9898

 visionflex.com.au



visionflex
Equality healthcare