

Accu-Chek Blood Glucose Meter

Wireless blood glucose monitoring system

Wireless connection to your Visionflex ProEX telehealth device or compatible smartphone, keeps your important information at your hand where and when you need it – eliminating the need for a paper diary. Along with proven accuracy, features like the spill resistant vial, wide test strip dosing area, strip port light and strip ejector button can help make every day glucose monitoring easier. The Accu-Chek Guide Blood Glucose Monitor has a measuring range of 0.6 mmol/L to 33.3 mmol/L to provide better management and tracking.

**Proven accuracy
that exceeds
current ISO
standard**



Accu-Chek Blood Glucose Meter

visionflex
Equality healthcare

Features & benefits

- **Proven accuracy that exceeds current ISO standard**
- **Connects wirelessly to Visionflex ProEX telehealth devices**
- **Have the option of transmitting results to the Accu-Chek Connect diabetes management app on your compatible smartphone using wireless technology**
- **No more “strip confetti” – the spill resistant vial helps you remove one strip at a time**
- **Wide test strip dosing area – apply a small drop of blood anywhere along the wide yellow edge**
- **Receive glucose readings in less than 4 seconds**
- **Monitor in all light conditions, thanks to a bright, backlit display and a strip port light**
- **Eject used test strips from the meter with the push of a button**
- **Enjoy virtually pain free lancing with the included Accu-Chek FastClix lancing device**

Specifications

Display – LCD with backlight

Blood sample size – 0.6µL

Power supply – Two 3-volt coin-cell batteries (CR2032)

Auto power-off after 90 seconds

Memory capacity – 720 blood glucose diary records plus 32 control tests

Dimensions – 80 x 47 x 20 mm

Bluetooth wireless technology

Data port – Micro USB (cable not supplied)

Test principle – Electrochemical

Test averages – 7, 14, 30, 90 days

Data management - On-board logbook, averages and pattern detection; Connects wirelessly to the Accu-Chek Connect diabetes management app; Data is stored in the cloud and can be accessed online via the Accu-Chek Connect Online system

System operating conditions – 2°C - 30°C (36°F - 84°F)

For further information and a demonstration, contact **Visionflex** at info@visionflex.com.au

1300 059 926 +61 2 8438 9898

visionflex.com.au



visionflex
Equality healthcare