

Flexible Video Rhino Laryngoscope

visionflex
Equality healthcare


3.4mm insertion tube diameter

High-definition flexible video rhino laryngoscope for visual examination of upper respiratory tract structures, including the nasal tract, nasopharynx, oropharynx and larynx.

Featuring built-in LEDs for enhanced illumination, and multi-function buttons. Offers 100° field of view with up to 140° bending angle.



HD images delivering detailed observations.



Flexible Video Rhino Laryngoscope

3.4mm insertion tube outer diameter

Flexible Video Rhino Laryngoscope with CMOS image sensor delivers high-definition quality images of nasal cavity, pharynx and larynx.

Featuring built-in LEDs for enhanced illumination, and multifunction buttons to control white balance, brightness and freeze functionality.

Offering 100° field of view with up to 140° bending angle, have confidence that you are getting the best, close-up view of your patient.

IN THE BOX

- Flexible video rhino laryngoscope
- USB video adapter cable
- Waterproof cap
- Portable air pressure leak detector
- User manual

FEATURES & BENEFITS

- Slim insertion tube measuring just 3.4mm
- CMOS image sensor delivers high quality imagery
- Built in LEDs provide illumination
- High-definition image enables more detailed observation of nasal cavity, pharynx and larynx
- Multifunction buttons to control white balance, brightness and freeze functionality
- USB output, UVC driver free
- 300mm insertion tube
- 100° field of view
- Up to 140° bending angle of the bending section

SPECIFICATIONS

Resolution	1280 x 720 HD
Insertion tube	3.4mm
Distal end	3.6mm
Working channel	No
Bending angle (up/down)	140°/140°
Working length	300mm
Connection	USB Connector, directly to PC or Laptop, UVC driver free, compatible to all kinds of image processing software
Operating conditions	Ambient temperature: 5°C - 40°C Relative humidity: 30% - 85%
Storage conditions	Ambient temperature: -20°C - 50°C Relative humidity 5% - 95%



For more information or to book a demonstration, please contact Visionflex at sales@visionflex.com



+61 2 8914 4000



visionflex.com



visionflex
Equality healthcare