

Horus Digital Scope HD

visionflex
Equality healthcare


Multi integrated HD imaging system with interchangeable lenses.

The Horus Digital Scope HD is an easy-to-use hand-held video system used for capturing videos and images of the body. This system utilises interchangeable attachments making it ideal for multi discipline medical applications.

Attachments for ear, nose and throat, dermatology, ophthalmology, audiology, women's health and general examinations are just a few of the potentials of this unique system.



Touch screen operation.



Horus Digital Scope HD

Multi integrated HD imaging system with interchangeable lenses.

The Horus Digital Scope HD incorporates high resolution camera technology and offers multiple video output options for connection to Visionflex ProEX telehealth devices or conferencing system. Still images and videos can be captured with the touch of a button and transferred seamlessly for storage and forwarding purposes.



FEATURES & BENEFITS

- Touch screen operation
- Manual focus with live image zoom
- Image capture button
- Brightness adjustment
- Interchangeable lens attachments are kept secure with the lens lock (Series 3)
- HDMI and USB live video connections with lockable cables (Series 3)



SPECIFICATIONS

Resolution

1920 x 1080p HD

Light Source

Light emitting diode (LED)

Display

3.5" / 90mm full colour LCD

Image Format

JPEG (still) & H.264 (video)

Media Card

Micro SD 2GB (included)

Power Source

Rechargeable Lithium Ion Battery

Operating Time (approx.)

3 hours

Weight

340g

Focusing Range

5mm to 30mm

SERIES 1

SERIES 3

2560 x 1920p HD (5Mp)

Light emitting diode (LED)

3.5" / 90mm full colour LCD

JPEG (still) & H.264 (video)

Micro SD 16GB (included)

Rechargeable Lithium Ion Battery

3 hours

275g

7mm to 20mm

For more information or to book a demonstration, please contact Visionflex at sales@visionflex.com

+61 2 8914 4000

visionflex.com



visionflex
Equality healthcare