


# ProEX Mobile



Portable, rugged, battery powered all-in-one telehealth tablet.

The ProEX Mobile is a portable, rugged, battery powered all-in-one telehealth tablet that integrates clinical examination functionality with video conferencing to deliver the next generation of remote medical care.

ProEX Mobile allows health practitioners to carry out detailed medical examinations on patients based in regional, remote, and urban areas; host live video consultations; connect with colleagues and specialists over large distances; and share saved images and clinical data.



**Connection to a specialist for remote examinations**



# ProEX Mobile

Portable, rugged, battery powered all-in-one telehealth tablet.

Vital patient data can be collected and viewed live on screen, and recorded and saved in the internal patient database.

Using VF Sync, multiple ProEX units can be easily programmed to communicate securely with each other, and information can be shared in a customised network location or uploaded to an electronic health record (EHR) system.

With store-and-forward capabilities, ProEX Mobile is perfect for fieldwork, outreach clinics and travelling health professionals in rural and remote locations including remote worksites. All patient data is encrypted, ready to comply with HIPAA and GDPR requirements.

Perform remote general imaging; ENT examinations; wound management; dental examinations; ophthalmology; dermatology; mental health; speech pathology and more. Other patient observations can be conducted using Bluetooth® and USB-connected devices.



## FEATURES & BENEFITS

- Rugged, state of the art tablet with kick stand
- Dual hot-swappable rechargeable battery design
- Super bright, sunlight readable 11.6" 16:9 LCD with touch screen
- IP54 water/dust rated shock proof housing
- Powerful patient data management software and archiving system
- High-quality live video and imaging capabilities
- HIPAA and GDPR ready
- Integration-ready with major electronic health record systems (EHRs)
- Integrated video conferencing software with webcams and microphones built-in front and back
- Conduct visual medical examinations such as general imaging, wound management, dental, ophthalmology, dermatology and more
- Perform clinical examinations with Bluetooth® devices such as ECG, pulse oximeter, blood pressure monitor, thermometer and blood glucose meter
- Optional: tabletop and vehicle docking stations

For more information or to book a demonstration, please contact Visionflex at [sales@visionflex.com](mailto:sales@visionflex.com)

+61 2 8914 4000

 [visionflex.com](http://visionflex.com)



**visionflex**  
Equality healthcare