

# Professional Extra-Large Weighing Scales

Comfortable, stable, and precise scales for weight control. Ideal for medical use.

These SECA extra-large scales have been designed in Germany to meet the highest demands in regard to comfort, stability and precision. SECA is the world market leader for medical weighing systems.

Connect to a Visionflex ProEX telehealth device or smart phone via Bluetooth® and the measured data will be transmitted to the phone synchronously.

Compatible with Google Fit and Apple Health.



Wireless connection  
via Bluetooth®.



# Professional Extra-Large Weighing Scales



Comfortable, stable, and precise scales for weight control. Ideal for medical use.

## FEATURES AND BENEFITS

- Wide footprint scales for all body shapes and sizes
- Robust steel construction with non-slip rubber feet and a soft, pimped surface for comfort
- 100g increments
- Large and clear LCD display
- Bluetooth® integration
- Auto power off
- German designed
- Requires 4 x AAA batteries (included)
- CE Certified

## SPECIFICATIONS

<b>Maximum weight capacity</b>	200kg (440lb)
<b>Measurement increments</b>	0.1kg
<b>Accuracy</b>	±0.2kg
<b>Measurement unit</b>	kg, lb or st
<b>Dimensions</b>	433 x 373 x 47mm
<b>Net weight</b>	2.9kg
<b>Power source</b>	4 x AAA batteries (included)
<b>Display</b>	LCD screen
<b>Operating temperature</b>	10°C - 40°C
<b>Bluetooth®</b>	4.0 technology



For more information or to book a demonstration, please contact Visionflex at [sales@visionflex.com](mailto:sales@visionflex.com)

 +61 2 8419 4000

 [visionflex.com](http://visionflex.com)



**visionflex**  
Equality healthcare